



We Help You To Be Safe



+20 Years Experience In The Medical Field

**Supporting Patient Safety And
Supply New And Unique Products**

20th Years With experience in the following



- ✓ Dialysis
- ✓ Cardiothoracic & Vascular surgery
- ✓ Infusion Therapy & TPN
- ✓ Anaesthesia and ICU
- ✓ Dermatology
- ✓ Clinical Nutrition
- ✓ Oncology
- ✓ Interventional Cardiology
- ✓ General surgery
- ✓ Critical care & operation
- ✓ Infection control
- ✓ Diabetic Care

Critical Care

The line is divided into 3 categories:

Critical Care

From innovative hemodynamic monitoring and fluid management to essential tools that help you better manage your patients' airways, temperature, and pain.

Oncology and infusion therapy

Help keep pharmacists, nurses, and patients safe with the needle-free closed system transfer devices, port, and drug compounding systems for handling hazardous IV medications—designed to be easy to use, safe, effective, and affordable.

Infusion Therapy (Ameco's disposable)

- Anesthesia and ventilation breathing system.
- Oxygen and Aerosol therapy system

Products:

- **Clave family**
- **Transpac**
- **DOM**
- **Pulmonary Artery Catheter**
- **Diana system**
- **AMECO disposables**

Clave family

Microclave

- Help Reduce Risk of CRBSIs
- Visualize Connector Flushing
- Use on all Catheters
- No Change in Technique



Silicone Seal and Internal Cannula Minimizes Point of Entry for Bacteria

Specifically designed to minimize contact between the connector's external surface and the internal fluid path upon luer activation, Clave needle free IV connectors minimize entry points for bacteria. Several studies have attributed this feature to a reduction in bacterial contaminants passed through the connector.

Split Septum

Clave's normally closed swabbable split septum design is a preferred feature for needlefree connectors.

Straight Fluid Path

Clave's straight fluid path allows for efficient clearing of medications, blood, and blood residual with low flush

Minimal Residual Volume

Clave's minimal residual volume allows for lower flush volumes.

Clear Housing

A clear housing lets you see whether you have completely flushed the connector after blood draws or administration.





Clave needle free connector.




Clave needle free connector with extension.



Clave needle free connector with Double extensions.



Clave needle free connector with Triple extensions.



Technical Specifications

Residual Volume	0.04 ml
Flow Rate at Gravity	185ml/min
Functional Activations	600
Blood Compatibility	Yes
MRI Compatibility	Yes

Drug Compatibility

Alcohol	Yes
Lipids	Yes
Chlorhexidine	Yes
Chemotherapy	Yes

CemoClave family

CemoClave

- Needlefree Closed System Transfer Device (CSTD).
- Maintain a needlefree closed system to help minimize exposure to hazardous drugs and comply with USP <800>.
- Intuitive system with CSTD bonded to help you improve workflow efficiency and prevent the risk of accidental HD exposure due to disconnect.

Bag Spike




Vial Axis



Hemodynamics

Transpac IV

- Transpac IV Disposable Pressure Transducers provide clinicians with real-time access to their patients' hemodynamic and cardiovascular status in surgical and intensive care settings.





Technical Specs

Operating Pressure	500 to 300 mm Hg
Overpressure Limits	400 to 5000 mm Hg
Sensitivity	5 microvolt/volt/mm Hg \pm 1% @ 6 HDC and 22°C
Zero Offset	< 25 mm Hg
Zero Drift	< 2 mm Hg in 8 Hours
Input Impedance	300-350 ohms
Output Impedance	300 ohms Nomina

DOM

Measure the invasive blood pressure



TORQUE-LINE

TD Thermodilution catheter (Pulmonary Artery Catheter)

Designed to provide accurate Pulmonary Artery Pressure (PAP) and Pulmonary Artery Occlusion Pressure (PAOP or wedge pressure) monitoring.


- PA Temperature Measurement.
- PAP and PAOP Monitoring.
- Venous and Mixed venous blood sampling.
- Central Venous pressure (CVP) Monitoring.
- IV drug administration.
- 4 lumen, 7 French sizes & 110cm.
- Heparin coated.
- Made in USA and FDA approved.



DIANA Hazardous Drug Compounding System

- Accuracy, Safety, Efficiency, Sterility.
- User-controlled automated compounding for maximum accuracy & safety.
- Closed system assures the safety of clinicians and the sterility of the mix.
- Reduces risk of repetitive stress injuries.
- Increases efficiencies and reduces drug waste.

The Diana System at a Glance



Ameco


Anesthesia Breathing Systems




- ✓ All colors available
- ✓ Lightweight corrugated tube
- ✓ No leakage
- ✓ Economical systems
- ✓ Connector and tube tightly fitted
- ✓ Luer cap options
- ✓ Without glow
- ✓ Comfortable use and flexible
- ✓ No kinking risk




Corrugated Basic Systems



**Corrugated Basic Systems
with extra limp**




**Corrugated Basic Systems
with extra limp, breathing
bag**




**Corrugated Basic Systems
with extra limp, filter,
breathing bag**


Ventilation Breathing Systems




Ventilation circuits with extra limb, single water trap



Ventilation circuits with extra limb, double water trap



Ventilation circuits with extra limb, double water trap, bacterial filter



Ventilation circuits with extra limb, double water trap, T-piece nebulizer



Resuscitators




**Ventilation circuits for
portable ventilator**



**Ventilation circuits for
portable ventilator,
with bacterial filter**



**Anesthesia Mask
Air Cushion**



**Ventilation circuits with extra
limp, double water trap,
sampling line**



**Flexible fixed elbow
catheter mount**



**Flexible double swivel
catheter mount**



**Flexible straight
catheter mount**



Bacterial - Viral Filter

Oxygen Therapy Systems



**Oxygen Mask
(Luxury Series)**




**Oxygen Mask
(Elegant Series)**



**Oxygen Mask
(Eco Series)**



**Non-Rebreathing
Mask**



Nasal Cannula



**Venturi Mask with
corrugated tube**



**Venturi Mask without
corrugated tube**



**Venturi Mask
Fixed volume 40%**



Recovery Kit 40%

Aerosol Therapy System



**Nebulizer Mask
(Luxury Series)**



**Nebulizer Swivel
(Elegant Series)**



**Nebulizer Mask
(ECO Series)**



Atomizer Chamber



Nebulizer Mouth Piece



Nebulizer Mouth Piece



Nebulizer T-Piece

Dialysis Line

We offer a wide range of AMECATH high-quality catheters for short-term, long-term, and peritoneal use, each designed with patient comfort, safety and ease of use in mind. In addition to a wide range of bloodlines and Av fistula.

Products:


- **Short term hemodialysis catheter**
- **Long term hemodialysis catheter**
- **Peritoneal catheter**
- **Fistula**
- **Blood lines**

Short Term Haemodialysis Catheter

An Optimum Compromise between Competing Factors.

Where a balance between:

- Flexibility & and firmness.
- Softness and anti-kinking.
- Material and durability.
- Size and flow.



Ordering Information

Fr. Size	Length
06	08Cm
07	10Cm
07	13Cm
08	10Cm
08	12Cm
08	13Cm
09	10Cm
09	12Cm
09	13Cm
09	13Cm
09	15Cm
10	10Cm
10	13Cm
10	15Cm
12	12Cm
12	13Cm
12	15Cm
12	16Cm
12	17Cm
12	20Cm
12	24Cm
14	15Cm
14	17Cm
14	20Cm
16	15Cm
16	17Cm
16	20Cm

Kit Contents

Standard Kits	Standard Kits with Guiding syringe
Tapered Tip Dual Lumen Haemodialysis	Tapered Tip Dual Lumen Haemodialysis
J End Nitinol Guide Wire with Scaled Dispenser	J End Nitinol Guide Wire with Scaled Dispenser
Vessel Dilators	Vessel Dilators
Syringe	Guiding Syringe
Introducer Needle, Echogenic	Introducer Needle, Echogenic
Scalpel	Scalpel
Injection Cap	Injection Cap
Movable Fixation Wings (with pre-curved catheter)	Movable Fixation Wings (with pre-curved catheter)
Guiding Y Connector	

Step Tip Double Lumens Catheters Kit (SDLC)



Step Tip Double Lumens Catheter

Step Tip Design

- For maximum flow
- Decrease Recirculation

Kit Contents

Standard Kits	Standard Kits with Guiding syringe
Step Tip Soft Dual Lumen Haemodialysis	Step Tip Soft Dual Lumen Haemodialysis
Stiffener with Haemostatic Valve	Stiffener with Haemostatic Valve
J End Nitinol Guide Wire with Scaled Dispenser	J End Nitinol Guide Wire with Scaled Dispenser
Vessel Dilators	Vessel Dilators
Vessel Dilator with hydrophilic coating	Vessel Dilator with hydrophilic coating
Syringe	Guiding Syringe
Guiding Y Connector	Introducer Needle, Echogenic
Introducer Needle, Echogenic	Scalpel
Scalpel	Injection Cap
Injection Cap	Unifix Catheter Tube Fixation Adhesive
Unifix Catheter Tube Fixation Adhesive	Transparent Catheter Fixation Adhesive
Transparent Catheter Fixation Adhesive	Movable Fixation Wings
Movable Fixation Wings	

Ordering Information

Fr. Size	Length
12	12Cm
12	15Cm
12	17Cm
12	20Cm
12	25Cm
14	15Cm
14	17Cm
14	20Cm
14	25Cm
14	30Cm
14	35Cm

Unicath Dual Lumens Catheters Kit (UN-DLC)

- Soft material: Comfortable for the patient & easily inserted.
- Smart groove: Protect side holes from sticking to vessel wall.
- Tubular design: No turbulence, higher flow.
- No connection: Prevent turbulence, more comfort.
- No clamp: Long term saved integrity of extensions.
- Hemostatic valves: Prevent blood leakage all over the procedure.



Kit Contents

Standard Kits	Standard Kits with Guiding syringe
Dual Lumen Haemodialysis Catheter	Dual Lumen Haemodialysis Catheter
Stiffener with Haemostatic Valve	Stiffener with Haemostatic Valve
J End Nitinol Guide Wire with Scaled Dispenser	J End Nitinol Guide Wire with Scaled Dispenser
Vessel Dilator	Vessel Dilator
Vessel Dilator with hydrophilic coating	Vessel Dilator with hydrophilic coating
Syringe	Guiding Syringe
Guiding Y Connector	Introducer Needle, Echogenic
Introducer Needle, Echogenic	Scalpel
Scalpel	Injection Cap
Injection Cap	Unifix Catheter Tube Fixation Adhesive
Unifix Catheter Tube Fixation Adhesive	Transparent Catheter Fixation Adhesive
Transparent Catheter Fixation Adhesive	Movable Fixation Wings
Movable Fixation Wings	

Ordering Information

Fr. Size	Length
12	15Cm
12	17Cm
12	20Cm
12	25Cm
14	15Cm
14	17Cm
14	20Cm
14	25Cm



Smart Grooves Design Prevents the side Holes from being Stick to the vessel wall During blood Suction.

Triple Lumens Catheters Kit (TLC)

Third lumen gives greater flexibility for fluid administration, vasoactive drug therapy, nutrition and blood sampling.



Kit Contents

Standard Kits	Standard Kits with Guiding syringe
Dual Lumen Haemodialysis Catheter	Dual Lumen Haemodialysis Catheter
Stiffener with Haemostatic Valve	Stiffener with Haemostatic Valve
J End Nitinol Guide Wire with Scaled Dispenser	J End Nitinol Guide Wire with Scaled Dispenser
Vessel Dilator	Vessel Dilator
Vessel Dilator with hydrophilic coating	Vessel Dilator with hydrophilic coating
Syringe	Guiding Syringe
Guiding Y Connector	Introducer Needle, Echogenic
Introducer Needle, Echogenic	Scalpel
Scalpel	Injection Cap
Injection Cap	Unifix Catheter Tube Fixation Adhesive
Unifix Catheter Tube Fixation Adhesive	Transparent Catheter Fixation Adhesive
Transparent Catheter Fixation Adhesive	Movable Fixation Wings
Movable Fixation Wings	

Ordering Information

Fr. Size	Length
12	15Cm
13	15Cm
13	17Cm
13	20Cm
13	24Cm
13	25Cm
14	15Cm
14	17Cm
14	20Cm

Long-Term Haemodialysis




Intended for implantable,
long-term central venous hemodialysis

Kit Contents

PLCD	
COMPONENT	INCLUDED
Long term dual lumen haemodialysis catheter	✓
J End Nitinol Core Guide Wire	✓
Stiffener with Cap	✓
Vessel Dilator	✓
10 cc Syringe	✓
G18 x 7 cm Introducer needle, echogenic	✓
Scalpel # 11	✓
Valved Peelable introducer sheath	✓
Blunt Curved Tunneling Stylet	✓
Injection cap	✓
Transparent catheter fixation adhesive	✓
Guiding Y-Connector	✓


Ordering Information

Fr. Size	Length
12	21Cm
12	24Cm
14	24Cm
14	28Cm
14	32Cm
14	36Cm
14	40Cm
14	55Cm
16	24Cm
16	28Cm
16	32Cm
16	36Cm
16	40Cm
16	55Cm




Step Tip
With Side
Holes

Peelable
introducer
sheath



Tunnelling Stylet



Cuff

"PYDLC" Configuration Permthane long-Term Y-Split Dual Lumen Haemodialysis Catheter Kit




Kit Contents

PLCD	
COMPONENT	INCLUDED
Long term Y-split dual lumen haemodialysis catheter	√
J End Nitinol Core Guide Wire	√
Stiffener with Cap	√
Vessel Dilator	√
10 cc Syringe	√
G18 x 7 cm Introducer needle, echogenic	√
Scalpel # 11	√
Valved Peelable introducer sheath (Peelable introducer sheath for 16Fr catheter)	√
Blunt Curved Tunneling Stylet	√
Injection cap	√
Transparent catheter fixation adhesive	√
Guiding Y Connector	√

Ordering Information

Fr. Size	Length
12	21Cm
12	24Cm
14	24Cm
14	28Cm
14	32Cm
14	36Cm
14	40Cm
14	55Cm
16	24Cm
16	28Cm
16	32Cm
16	36Cm
16	40Cm
16	55Cm

“PXDLC” Configuration Permthane X-Split long-Term Dual Lumen Central Venous Hemodialysis Catheter Kit



Kit Contents

PLCD	
COMPONENT	INCLUDED
Long-Term X-split dual lumen haemodialysis catheter	√
J End Nitinol Core Guide Wire	√
Movable fixation wing	√
Vessel Dilator	√
10 cc Syringe	√
G18 x 7 cm Introducer needle, echogenic	√
Scalpel # 11	√
Valved Peelable introducer sheath (Peelable introducer sheath for 16Fr catheter)	√
Blunt Curved Tunneling Stylet	√
Free female luer line for retrograde (Red for arterial + Blue for venous)	√
Flushing connector	√
Slide clamps	√
Catheter tube fixation adhesive (Unifix)	√
Transparent catheter fixation adhesive	√
Guiding Y-connector	√

Ordering Information

Fr. Size	Length
12	21Cm
12	24Cm
14	24Cm
14	28Cm
14	32Cm
14	36Cm
14	40Cm
14	55Cm
16	24Cm
16	28Cm
16	32Cm
16	36Cm
16	40Cm
16	55Cm

Long-Term Peritoneal Dialysis Catheter Kits



Kit Contents

PSP (Straight Peritoneal Catheter Kit)	
COMPONENT	INCLUDED
Straight peritoneal dialysis catheter with one or two cuffs	√
J-End Nitinol core 0.035" x 70 cm guide wire	√
Vessel Dilator	√
10 cc Syringe	√
G18 x 7 cm Introducer needle, echogenic	√
Scalpel # 11	√
Grip-Lock peelable introducer sheath	√
Blunt curved tunnelling stylet 18cm	√
Movable fixation wing	√
Haemostatic valve	√
Transparent catheter fixation adhesive	√
Guiding Y-connector	√
Free female luer line for retrograde-blue (for cath. size 10Fr only)	√
Female luer (for cath. size 12, 14 and 16Fr)	√
Clamp (for cath. size 12, 14 and 16Fr)	√

Ordering Information

Fr. Size	Length
10	1Cm
10	1Cm
10	2Cm
12	1Cm
12	2Cm
12	2Cm
14	1Cm
14	1Cm
14	1Cm
14	1Cm
14	2Cm
14	2Cm
16	1Cm
16	1Cm
16	2Cm
16	2Cm

Pediatric Long-term hemodialysis catheter



Kit Contents

Pediatric Long-term hemodialysis catheter	
COMPONENT	INCLUDED
Long-Term Pediatric Twin lumen Haemodialysis catheter kit with integrated Female Luer Valve	√
J-End Nitinol core guide wire	√
Vessel Dilator	√
5 cc Syringe	√
G18 x 5 cm Introducer needle, echogenic	√
Scalpel # 11	√
Movable fixation wing	√
Valved peelable introducer sheath	√
Blunt curved tunneling stylet	√
Repair valves	√
Transparent catheter fixation adhesive	√
Catheter tube fixation adhesive (Unifix)	√
Guiding Y-connector	√

Ordering Information

Fr. Size	Length
8	12Cm
8	15Cm
8	18Cm
10	12Cm
10	15Cm
10	18Cm
10	23Cm
12	28Cm
12	32Cm

AV Fistula Needles




Product Features


- Smooth tubing in-wall avoid blood coagulation and bubble.
- Anti-fatigue pump tube keeps same blood flow in continuous 50 hours.
- Assured product safety by permanent in-process controls and 100% leakage tests.
- A wide range of different configurations is available to suit all dialysis machines.
- Provided with self sealing injection sites for intermittent medication.
- Drip chamber of venous line contains filter to prevent clot entering into patient's vein.
- Transparent, soft and kink resistant tubing made of medical grade PVC.

Product variations

Component/ Aspect	Options
Cannula bevel type	Sharp multiangle Blunt
Cannula shaft (inch)	1"
Cannula backeye	With Without
Needle size (guage)	18G, 17G, 16G, 15G, 14G
Tube length (cm)	20, 25, 30, 35, 40
On/Off clamp	White Blue Red
Safety guard	With Without
Packaging	Single Double



Blood Lines



Product Features

- Smooth tubing in-wall avoid blood coagulation and bubble.
- Anti-fatigue pump tube keeps same blood flow in continuous 50 hours.
- Assured product safety by permanent in-process controls and 100% leakage tests.
- A wide range of different configurations is available to suit all dialysis machines.
- Provided with self sealing injection sites for intermittent medication.
- Drip chamber of venous line contains filter to prevent clot entering into patient's vein.
- Transparent, soft and kink resistant tubing made of medical grade PVC.

Product variations

Component/ Aspect	Options
Main tubing diameter	6.8*4.6mm
Tubing length venous side (cm)	275, 300
Tubing length arterial side	270, 300
Pump segment diameter	12*8mm 9.8*6mm
Pump segment length	33cm
Drip chamber type	22/30
Drip chamber cover	2-way 3-way 4-way
Drip chamber cover tubing	3.5*5.5mm
Injection site	2-way 3-way
Standard components	Dialyzer connector & cap Female luer & lock ON/OFF bench clamp Pressure transducer Rotating male luer & lock
Optional components	Priming bag Recirculating Hook
Device type compatability	Fresenius Gambro Bbraun Nekisso Baxter



Infection Control

Safehearts Provide integrated solution for Disinfection, Sanitization, and Sterilization for a healthy environment in Health Care Facilities and public Environments.

We provide a Complete range of disinfection solutions for treatment:

- Surfaces, Hands, Tools, instruments, Patient skin, and Environment (Air and surfaces)

Products:

- Novaerus Air Sanitization Systems**
- Rein products for healthy Environments**
- NADES Chemical disinfectants**
- Hexa Chemical disinfectants**

NOVAERUS

Novaerus Cleans The Air & stops the infection circle with the NanoStrike technology


Advantages & benefits

- The Only devices FDA, 510K, and CE Approved (Insured quality, safety, and effectiveness).
- Nano-time killing of all airborne Pathogens with Nano strike technology.
- Best CADR (Clean Air Delivery Rate) in the Market Up to 96.25 %.
- Achieving clean room standard Up to class 7 clean room ISO Standard for a healthy environment.
- Achieve Best Room ACH (Air Chang /Hour) up to 4 ACH.
- Protection up to 99.97% (HEPA Filter Camfile H13), According to NOISH STANDARD.
- Complete removal of all chemicals, odors & and VOCs by Carbon Filter.

Multiple designs for different requirements

NV protect 200

- Dimensions: 28.3 cm (h) x 13.2cm (w) x 10.8 (d).
- Electric consumption: 14 W
- 1 NanoStrike coil.
- WEIGHT 3.4kg | 7.5lbs



NV protect 800

- Dimensions: 36.5 cm (h) x 36.5 cm (w) x 11.5 cm (d).
- Electric consumption: 25 W
- 2 NanoStrike coils.
- FILTERS
- Stage 1 Pre Filter
- WEIGHT 4.7kg | 104 lbs.



NV Defend 400

- Dimensions: 85.5 cm (h) x 38.5 cm (w) x 38.5 cm (d).
- 4-NanoStrike coils,
- Stage 1 Pre Filter: M5 EN 779, ePM10 65% ISO 16890
- Stage 2 HEPA Filter: H13 EN 1822
- Stage 3 Carbon Filter: G4 EN 779, Coarse 65% ISO 16890
- WEIGHT 25.7kg | 56.7 lbs.



NV Defend 1050

- Dimensions: 93 cm (h) x 58 cm (w) x 48 cm (d).
- POWER CONSUMPTION RANGE
- Lowest Speed: 60.5W – 78.5W
- Highest Speed: 309.7W – 395.8W
- 6-NanoStrike coils,
- 1 Pre Filter: M5 30/30 filter
- Stage 2 HEPA Filter: HEPA H13 filter, MGA grade
- Stage 3 Carbon Filter: G4 carbon pleated filter
- WEIGHT 53.8kg | 118.6 lbs.



Miller

Miller disinfectants are an excellent choice for disinfecting both contents and surfaces with far lower concentrations.

Miller Disinfectants Are Effective Against

- Gram-positive and Gram-negative bacteria.
- viruses, including bacteriophages, spore-forming microorganisms, yeast, fungi, mold & protozoa.

- ✓ Highly effective biofilm remover
- ✓ Universal stability
- ✓ Forms no harmful by-products
- ✓ Non-toxic and Ph-neutral when dosed
- ✓ Proven efficiency and globally approved
- ✓ Odor, color, and taste-free
- ✓ Cover all hidden areas in the room
- ✓ Very short contact time (30-60 min)
- ✓ Cost effective (utilize very small volume)
- ✓ Elimination of any human factor



Pack Volume: gallon 5 liter
Concentration: 3%, 6%, and 12%

Safe Fog Rotator

Electric **ULV** Cold Fogger



Why Safe Fog Rotator is the Ideal fogger for the safest and most effective disinfection solution **miingst** in the big rooms up to 3000 m³

Advantages

- Highly versatile and ideal equipment for the application of disinfectants in all health care facilities, hospital and also in public places.
- Range of the fog is deeper than with comparable units.
- Range of application includes space and surface treatment
- Easy handling and operating of the unit

Rein detergent

Anti-Bacterial Hand Soap

- ① Reduce Load of contamination and organic materials.
- ② Recommended for hygienic hand wash in food industry.
- ③ Ready to use.
- ④ Leave hand fresh & clean.
- ⑤ It cleans dirt, oils, grease and other unwanted substances from hand.
- ⑥ Cleaning and pre-sanitation in one step.
- ⑦ For hygienic hand wash in hospitals, medical sector, food industry, and pharmaceutical industry.
- ⑧ Liquid.
- ⑨ Odorless.



No Smell



No Color

Detergent

Multi-Purpose Detergent Anti-bacterial

- ① Concentrated highly alkaline formula.
- ② Multipurpose low foam detergent.
- ③ Cleaning and pre-sanitation in one step.
- ④ Effective in hard and soft water.
- ⑤ Used in food, beverage and pharmaceutical industries.



Anti-Bacterial Chloro- Foam

Benefits:

- ① Specially developed for use with low pressure foam devices.
- ② Working solutions from stable foam that adheres for a long time, even to vertical walls and surfaces. Areas those are difficult to clean manually can thus be cleaned reliably and safely.
- ③ Cleaning and pre-satination in one step, where it reduces load of contamination, organic materials and helps to prevent build-up of protein film.
- ④ Assist stain removal.
- ⑤ Suitable for use on most types of surfaces, but unsuitable for use on soft metals and galvanized materials.
- ⑥ Effective in hard and soft water,



5 - 10 - 20 kg

NADES

NADES & NADES 2.0 USES

For Disinfection of all Surfaces, tools, scopes, instruments, medical devices, and Equipment.

Super oxidized water

Electrolyzing saline

New Generation of Antiseptic



NADES⁺
Medical Tools
Disinfectant



NADES 2.0⁺
Surgical Instruments
and Tools Disinfectant



NADES⁺
Medical Surface
Disinfectant



NADES⁺
Medical Equipment
Disinfectant

Hexa family

Hexa hand



HEXA ETHANOL 80% LIQUID

- Description

Alcohol-based for hygienic and surgical hand antiseptic.

Hand antiseptics used in Healthcare, High Concentrate from Ethyl Alcohol to be fast action for reducing bacteria and viruses that potentially can cause disease.

- Composition

Active ingredients Each 100gms contains:

- Ethyl Alcohol 80% v/v.
- Glycerol 1.45% v/v.
- Purified water USP.



HEXA ISOPROPYL 75%

- Description

Alcohol-based for hygienic and surgical

Hand antiseptics used in Healthcare, High Concentrate from Alcohol to be fast action for reducing bacteria and virus that potentially can cause disease.

HEXA ISOPROPYL 75% is an Alcohol-based ready-to-use solution for hygienic and surgical hands.

- Composition

Active ingredients Each 100gms contains:

- Isopropyl Alcohol 75% v/v.
- Glycerol 1.45% v/v
- purified water USP.



HEXA 70% LIQUID

- Description

Alcohol-based for hygienic and surgical hand antiseptic
ETHYL ALCOHOL AND CHLORHEXIDINE GLUCONATE

- Composition

Active ingredients Each 100gms contains:

- Ethyl Alcohol 70%
- Chlorhexidine gluconate less than 1%
- Glycerin 3%
- Triethanolamine 3%
- purified water.



HEXA 70% GEL

- Description

Alcohol-based for hygienic and surgical hand antiseptic
ETHYL ALCOHOL AND CHLORHEXIDINE GLUCONATE

- Composition

Active ingredients Each 100gms contains:

- Ethyl Alcohol 70%
- Chlorhexidine gluconate less than 1%
- Glycerin 3%
- Triethanolamine 3%
- purified water.

Instruments and Equipment Disinfection Products



HEXA RUST PLUS

- Description

- Rust and Stain remover for surgical instruments, carts and trays
- Mild acidic solution removes rust, stains, spotting and corrosion
- Restores instruments original finish Non-corrosive; will not harm stainless steel
- Save money by restoring instruments instead of replacing or repairing them
- Routine use.

- Composition

Each 100gms contains:

- 30 g phosphoric acid
- <5 g nonionic surfactant
- Corrosion inhibitor
- Auxiliary substances.



HEXAZYM 5

- Description

Concentrated, 5 Enzyme-Based Liquid Cleaner
Economical and Fast-Acting Enzymatic Detergent

- Composition

- Non-ionic and amphoteric surfactants
- Solvents
- Complexing agents
- Amino acid derivatives
- Multi Enzyme complex (Amylase, Cellulose, Lipase, Mannanase and Protease).
- Corrosion inhibitors.
- Preservative.

SURFACES DISINFECTANT



HEXA SURF PLUS

- Description

- Detergent and disinfectant concentrated For Surfaces and Medical Devices
- Cleaning and disinfection for surfaces and Medical Devices such as Floors, walls, incubators, Equipment, etc.

- Composition

- 25 gm/g Didecylidimethylammonium chloride
- 50 gm/g N-(3-aminopropyl)-Ndodecylpropane1,3-diamine
- 6 gm/g Diethylene glycol monobutyl ether, perfume and colouring.
- Deionized Water.



HEXA SURF FOAM

- Description

- Detergent & Disinfectant Foam for Surfaces and Medical Devices

- Composition

- 0.25 Alkyl Ethylbenzyl Dimethyl Ammonium Chloride
- 0.25 Didecylidimethylammonium Chloride
- 0.25 Alkylbenzyl dimethyl Ammonium Chloride
- Auxiliary Substances and Deionized Water.



HEXA SURF RS

- Description

Alcohol-Based Rapid Spray Detergent disinfectant for surfaces and Medical Devices

- Composition

- 55% Ethanol
- 0,05% Didecylidimethylammonium chloride
- Alkaline agents Glycol derivatives
- Auxiliary Substances
- Deionized Water.

SURFACES DISINFECTANT



HEXA SURF PLUS

- Description

- Detergent and disinfectant concentrated For Surfaces and Medical Devices
- Cleaning and disinfection for surfaces and Medical Devices such as Floors, walls, incubators, Equipment, etc.

- Composition

- 25 gm/g Didecyldimethylammonium chloride
- 50 gm/g N-(3-aminopropyl)-Ndodecylpropane1,3-diamine
- 6 gm/g Diethylene glycol monobutyl ether, perfume and colouring.
- Deionized Water.



HEXA SURF FOAM

- Description

- Detergent & Disinfectant Foam for Surfaces and Medical Devices

- Composition

- 0.25 Alkyl Ethylbenzyl Dimethyl Ammonium Chloride
- 0.25 Didecyldimethylammonium Chloride
- 0.25 Alkylbenzyl dimethyl Ammonium Chloride
- Auxiliary Substances and Deionized Water.



HEXA SURF RS

- Description

Alcohol-Based Rapid Spray Detergent disinfectant for surfaces and Medical Devices

- Composition

- 55% Ethanol
- 0,05% Didecyldimethylammonium chloride
- Alkaline agents Glycol derivatives
- Auxiliary Substances
- Deionized Water.

Diabetic Care

We provide Sinocare biosensor technology, revolutionizing diabetic diagnosis with its innovative capabilities.

Products:

- **Safe AQ Max I**
- **Safe AQ Max III**
- **Safe AQ UG**
- **iCan**



Sinocare
◆ caring for love

Safe AQ Max I

Blood Glucose Monitoring System



Specification

FAD-GDH system

Blood Sample: fresh capillary blood, venous blood

Calibration Sample: venous blood

Test Time: 5s

Blood Volume: 1 μ L

HCT: 10%~70%

Test Range: 20~600 mg/dL (1.1~33.3 mmol/L)

Expiry Date of Test Strip: 24 months

Test Condition: 10°C~40°C

Meter Storage Temperature: -20°C~55°C, \leq 95%RH

Battery Life: above 1000 tests

Power: 2 AAA alkaline batteries 3V DC

Memory: 500 blood glucose test results,
100 control solution test results

Safe AQ Max III

Blood Glucose Monitoring System



Specification

FAD-GDH System

Blood Sample: Fresh capillary blood, venous blood

Calibration Sample: Venous blood

Test Time: 5s

Blood Volume: 1 μ L

HCT: 10%~70%

Test Range: 20~600 mg/dL (1.1~33.3 mmol/L)

Expiry Date of Test Strip: 24 months

Test Condition: 10°C~40°C

Meter Storage Temperature: -20°C~55°C, \leq 95%RH

Battery Life: above 1000 tests

Power: 2 AAA alkaline batteries 3V DC

Memory: 500 blood glucose test results

Safe AQ UG

Blood Glucose & Uric Acid Monitoring System



Specification

Test Function	Blood Glucose	Uric Acid
Test Metho	FAD-GDH	Oxidized Ferrocence Derivative
Sample Type	Capillary whole blood, venous whole blood	Capillary whole blood, venous whole blood
Sample Volume	0.6µl	3µl
Test Range	20-600 mg/dL(1.1-33.3 mmol/L)	3.0-20.0 mg/dL(181-1188 µmol/L)
Test Time	5s	25s
Memory	500 test results	100 test results
Expiry Date Of Test Strip	24 months	18 months



iCAN

The world's first and only 15-day continuous monitoring system

iCan CGM offers you longer glucose monitoring for up to 15 days. Manage your diabetes with more confidence and less effort

Easier diabetes management

- No scanning.
- No fingerpicks.
- View your glucose readings and trends anytime and anywhere.

Customizable glucose alerts

Set low and high glucose alerts that suit your target glucose level and keep you informed of hyper and hypoglycemia.

More insights into your glucose readings

Offer multiple glucose analysis to easily know your glucose trends and fluctuation.

- Today
- Multiple days
- Daily average
- Comparison

Adaptable to various scenes of life

Incredibly compact and weight only 6g.

Water resistant to 8 feet for up to 2 hours.

You can exercise, shower, swim, or even surf at ease.

Derma line

Our diverse range of Nano-static products addresses your skin and hair concerns with precision and efficiency

Nanostetic Mesotherapy:

- **Hair Capil** Vial
- **Rewnewnax** Vial
- **Melaglow** Vial
- **Dermaglow** Vial
- **Melasnow** Vial
- **Anadrox** Vial
- **Biosoft** Vial
- **Aquatone** Ampoule

Nanostetic Serums:

- **Top Aqua** Serum
- **Youth-Tone** Serum
- **Neutrax** Serum
- **Melashine** Serum
- **Melashine Plus** Serum
- **Melashine Extra** Serum

Xai Producats:

- **Anagenic-Pro**
- **Esscarox**
- **Intimate Care**
- **Sebio-Refiner**

Mesotherapy products

Hair Capil vial

The superior Hair formula



- Hair Capil is a Superior Formula fortified by the power of DHQG & EGCG

DHQG Taxifolin **EGCG-Glucoside**

- Improve hair Growth
- Moisturizer effect
- Nourishment effect



Renewmax vial

Botox like effect



By its antagonist effect on muscular nicotinic acetylcholine receptors which decrease muscle contraction.

- Skin Firming
- Enhance Extracellular Matrix Formation
- Nourishment & moisturizing



Melaglow vial

For deep **face** hyperpigmentation



- Decreases pigmentation and Whitening
- Antioxidant
- Moisturizer



Dermaglow vial

For treatment hyperpigmentation in sensitive areas



- Whitening of the sensitive Area.
- Exfoliator and moisturizing
- Antioxidant



Melasnow vial

Covers all Eye's Problems



- Reducing dark hallows by reducing melanin production
- Puffiness and Sagging reduction
- Stimulating synthesis of collagen and collagen fiber regeneration of the dermal matrix occur



Anadrox & Biosoft vial

Active Acne & Post Acne Scars



- Reduction of sebum secretion
- Exfoliator
- Anti-microbial action
- Anti-inflammatory action
- Moisturizer



Biosoft vial

A mixture of potent active ingredients has a superior role in the treatment of post-acne Scar



- Epithelialization & Anti-Fibrosis
- Enhancement of Extracellular matrix
- Exfoliator
- Astringent effect
- Healing properties and anti-oxidants
- Moisturizer

Aquatone Ampoule

Experience the Joy of Rejuvenation



- Skin Rejuvenation
- Astringent to decrease dilated pores
- Acts as a moisturizer



Serums

TOP AQUA Serum

Premium formula to provide a filler-like effect



- Palmitoyl tripeptide-38
- Centella Asiatica Extract
- Coenzyme Q 10
- Ascorbic Acid (Vitamin C)
- Collagen (marine source)
- Elastin
- Hyaluronic Acid



YOUTH-TONE Serum

Superior solution to treat aging marks



- Argireline
- SYN-AKE
- Collagen (marine source)
- Elastin
- Hyaluronic Acid
- Tocopherol
- Calendula Officinalis
- Folic Acid



Neutrax Serum

Ultimate soothing and repairing moisturizer



- Hyaluronic acid
- Caviar Extract
- Snail Extract
- Aloe Vera
- Apple Stem Cell
- Coenzyme Q-10



Melashine Serum

Ingredients of melashine peeling serums:

- Salicylic Acid
- Lactic Acid
- Malic Acid
- Glycolic Acid
- Mandelic Acid
- Aloe Vera
- Hyaluronic acid
- Panthenol

Melashine Concentrations

	Melashine	Plus	Extra
Skin Type	Sensitive and dry skin	Normal and mixed skin	Oily skin
Glycolic Acid	2.5%	3%	4%
Mandelic Acid	2.5%	3%	5%
Salicylic Acid	2%	2%	2%



Xai Products

Anagenic-Pro

For recurrent and established Chronic Hair Loss and thinning Follicles

6 MONTHS PROGRAMME

- Ampoules 5ml x 11Amp

(Intensive treatment phase). Once a day.
Every other day for 6 weeks

- Shampoo 200 ml

3 Times a week During the treatment period.

- Lotion 120ml

(6 : 8 puffs) once every other day on days ampoules are to be skipped for 6 weeks.
Then once everyday for 4.5 months.

Aminexil - Procapil® - Trichogen™ VEG
Saw Palmetto - Caffeine



Esscarox

Gel - 40 ml

- Use twice Daily for 6 months.
- Apply a fair amount on the scar area and spread the gel with your finger tips to create a uniform silicone layer over the scar.

Polydimethylsiloxane - Cyclopentasiloxane
Silicone Dioxide - Vitamin E



Intimate Care Lightening Cream

- Use twice daily for 2 months.
- And 5 days after laser treatment.

Alpha Arbutin - Lumiskin™ - Niacinamide
Kojic acid - Liquorice Extract
Hyaluronic acid - Lemongrass Extract



Sebio-Refiner

START AND END THE DAY THE RIGHT WAY

- Gel - 60 ml

Apply twice a day a small amount on clean facial skin covering all the mild to moderate acne flairs and black heads.

- Foaming Clenser - 150 ml

Use twice a day morning and evening
Apply a fair amount on the face
message the skin with your fingers
gently and then wash it away.

Salicylic acid 2% - Tea Tree oil 5%
Niacinamide 2% - Aloe vera Extract
Zinc Oxide



Clinical Nutrition

Your source for high-quality nutritional supplements designed to support your well-being and enhance your health. With a focus on premium ingredients and expert formulations, Bonpharma is committed to providing you with the tools you need to nourish your body and thrive.

Products:

- Avepan
- Dhaepal
- Flaka
- No Ros
- Spark



AVEPAN

What's AVEPAN:

Anti-stress Biological, Clinical Nano Nutrition for chemotherapy radiotherapy patients.

Indications:

Special clinical nutrition formula therapeutic, Biological, Nano nutrition For Cancer patients to improve general health conditions and support chemotherapy and Radiotherapy treatment.

Composition:

- 1- Fermented wheat germ extract & Nano selenium with micronized pectin carried on rich base antioxidants & biological active nutrient's.
- 2- 100% natural food.
- 3- Pineapple flavor .
- 4- Contain 98 nutrients.
- 5- Contain CO Q 10 as potent anti-inflammatory.

For:

Chemotherapy & Radiotherapy patients.

Dose:

1,2 tablespoons in the morning with 2 glasses (plenty) of water, juice , milk & yoghurt at once daily for maximum three months.

Storage:

Keep refrigerator

Packet:

500 g

**One packet AVEPAN will be:
used for min 20 days**

Research: Analog to Avemar.



DHEAPEL

What's DHEAPEL:

The best immune modulator that readjust, restore immune system function's and natural anti viral function throw complete nutrition formula & also 100% natural food.

Indications:

Special supportive therapeutic Nano nutrition For immunocompromised patients to improve general health conditions and support treatment.

Composition:

- 1- Natural Wheat germ with source of protein, yeast & lemon flavor.
- 2- Special formula of amino acids with vitamins, minerals, active Co-enzymes, redox factors and anti-oxidants mixed with wheat germ oil.
- 3-100 % Natural food.
- 4- lemon & mint & yeast flavor.
- 5- Contain 98 nutrients.
- 6- FDA Approved.

For:

ICU patient, Inpatient, Viral infection patients, Chemo and Radiotherapy patients.

Dose:

One or two tablespoons in 2 glasses (plenty) of water, juice or yoghurt.

Storage:

Keep refrigerator

Packet:

500 g

One packet from DHAEPAL will be used for min 20 days



NO ROS

What's NO ROS:

- Complete nutritional formula for ICU patients, critical ill patient and hepatic patient.
- Powerful antioxidant, anti stress & immune stimulant.

Indications:

- 1- ICU patients.
- 2- Inpatients inside the hospitals.
- 3- Hepatic & IBS patients.
- 4- Diabetic Patients.

Composition:

- 1- Natural wheat germ as a source of protein with pomegranate, sucralose & citric.
- 2- 100% natural food.
- 3- Contain 98 nutrients.

For:

ICU patient, Inpatient, Malnourished, critical ill patients, Hepatic Patients.

Dose:

1,2 tablespoons) in the morning with 2 glasses (plenty) of water, juice, milk & yoghurt at once daily for maximum three months.

Storage:

Keep refrigerator

Packet:

500 g

One packet No ROS will be used for min 20 days





SAFEHEARTS

We Help You To Be Safe



**2A New Naser City towers - Hassan Maa"moun
Extension St.Al Wafaa Wa-AlAmal Naser City - Cairo.**



www.safehearts.co



info@safehearts.co



+201023336789



+2220797177